

Durodent™

MOULD CHART
CARTA DE MOLDES



ACRYLIC TEETH
DIENTES ACRILICOS



Upper and lower anteriors
Superiores e inferiores
anteriores Mould No.



Upper and lower anteriors
Superiores e inferiores
anteriores Mould No.



Upper and lower anteriors
Superiores e inferiores
anteriores Mould No.

2D 7.7 9.5 42.0



4.6 8.7 31.5

3D 8.0 9.3 44.0



4.8 8.7 34.0

1N 8.2 9.9 44.5



5.0 8.9 34.0

2N 7.4 10.0 42.0



5.0 9.3 32.1

3N 8.0 10.0 44.5



5.0 9.3 34.1

3M 7.0 8.8 38.3



4.6 8.4 31.2

Baby



2P 8.5 10.9 46.5



4.9 9.8 34.5

3P 8.8 10.5 47.6



5.8 10.0 37.6

1H 8.5 11.8 47.5



5.7 10.6 37.8

4H 8.6 10.3 46.9



5.7 9.1 37.2

222 6.8 7.5 36.0



4.5 8.4 32.3

223 7.0 7.9 38.5



4.4 8.4 32.9

224 8.3 8.5 42.0



4.6 9.0 32.0

A25 8.8 9.5 46.5



5.2 9.3 35.0

A26 8.8 9.5 46.5



6.0 8.8 36.6

264 8.3 10.5 45.2



5.8 9.7 38.3

136 9.0 10.1 48.7



5.3 9.1 36.8

137 9.2 10.3 51.0



5.6 9.5 37.3

266 8.8 11.7 48.6



6.0 10.7 37.8

3R 10.1 12.6 54.0





Posteriors 33°
Posteriores 33°

Mould No.

30M

8.9 | 8.1 | 29.7



9.3 | 8.9 | 31.8

30L

8.9 | 8.6 | 29.0



9.7 | 9.0 | 30.9

32M

9.8 | 7.6 | 32.0



11.0 | 8.5 | 34.0

32L

9.5 | 8.0 | 31.0



10.6 | 8.6 | 34.7

34M

10.1 | 9.0 | 33.4




10.9 | 9.4 | 35.8

34L

10.1 | 9.3 | 33.4



11.1 | 9.4 | 34.3



Posteriors 20°
Posteriores 20°

Mould No.

29M

9.0 | 7.0 | 29.5



9.9 | 8.0 | 31.9

29L

8.9 | 8.0 | 30.0



10.0 | 8.0 | 31.5

31M

9.7 | 8.1 | 31.0



10.6 | 8.8 | 34.0

31L

9.3 | 8.4 | 30.9



10.5 | 9.0 | 34.4

32F

9.3 | 8.6 | 31.5



10.1 | 9.4 | 33.0



Posteriors 0°
Posteriores 0°

Mould No.

29Z

9.0 | 7.7 | 29.7



10.3 | 7.7 | 31.3

31Z

9.5 | 7.5 | 31.0



10.7 | 7.6 | 31.8

33Z

10.5 | 8.4 | 34.3



11.4 | 8.9 | 34.8



Functional
Funcional

Mould No.

Acrylic Teeth

Dientes Acrilicos Made in Medellin - Colombia

Anteriors (Anteriores)**Upper Superiores****Lower Inferiores**

2D	2D
3D	3D
1N	1N
2N	2N
3N	3N
3M	3M
2P	2P
3P	3P
1H	1H
4H	4H
222	222
223	223
224	224
A25	A25
A26	A26
264	264
136	136
137	137
266	266
3R	3R

Functional (Funcional)

32F

Posteriors (Posteriores)**Anatomical
33°****Semianatomical
20°****Flat plane
0°**

30L	29L	29Z
30L	31M	31Z
32L	31L	31Z
32L	29L	29Z
32L	31M	31Z
30M	29M	29Z
32L	31L	31Z
34M	31M	33Z
34L	31L	33Z
32M	31M	33Z
30M	29M	29Z
30M	29M	29Z
30L	29M	29Z
32M	29L	31Z
34M	31L	31Z
34M	31L	31Z
34L	29L	33Z
34L	31L	33Z
34L	31L	33Z
34L	31L	33Z

Chart of articulations**and combinations**Carta de Articulaciones
y combinaciones**Coronas****Referencias**

3D

3P

A25

A26

137

**Shade**

58

Shades (Colores)

Cross linked fluorescent
and traslucent